

Spernu

FRILARD

Procedure for Freezing Human Semen

The fresh sample is extended with four times its volume of a solution prepared by homogenizing fresh egg whites in a buffer of the following composition.

NaCl : 5.00 g/l
KCl : 0.20 g/l
CaCl₂ : 0.20 g/l
MgCl₂ · 6H₂O : 0.10 g/l
NaH₂PO₄ · H₂O : 4.00 g/l
NaHCO₃ : 2.00 g/l
glucose : 1.00 g/l
ion-exchange H₂O : to make 1000 ml.
pH⁺ : 6.5

The extender is made by adding this buffer to 25 cc fresh egg white until the vol. = 100 ml.
The pH⁺ of the extender is 7.4 - 7.6

1/2 hour } The extended ^{semen} sample is then slowly diluted with a solution containing 20g glycerol/100 ml. The solution is added stepwise ($\frac{1}{3}$ vol. of extended semen/6 min) until the volume is doubled. The final concentration of glycerol = 10g/100ml extended semen.

The glycerol carrier is made by adding a buffer of the following composition to 20 gm glycerol until the volume = 100 ml.

NaCl : 5.030 g
KCl : 0.20 g
CaCl₂ : 0.20 g
MgCl₂·6H₂O : 0.10 g
NaH₂PO₄·H₂O : 0.5 g
NaHCO₃ : 5.00 g
glucose : 1.00 g
H₂O to make 1000 ml

(p. H. 7.4)

The extended ^{sample} sample (sample with 10% glycerol) is then dispensed into ampoules. Ordinary Wasserman tubes 11-13 mm I.D. have been used. An ampoule made by plugging a 15 mm I.D. test tube with a Wasserman tube to give an annular space of 7.5-2.5 mm. affords more even cooling. The ampoules are sealed in the oxygen-gas flame, taking care that the sample is not heated. Except for the brief operation of sealing the samples are protected from rapid temperature fluctuations by immersion in a 100m temp 25-28°C water bath.

An equilibration time of 30 min is allowed after the last glycerol addition. Then the ampoules are immersed in a freezing mixture of the proper temperature. The baths used were the following:

-160°C	=	Liquid-solid isopentane Eutectic mixture of EtOH-H₂O
-80°C	=	CO ₂ - EtOH ^{in mixture} 73% CO ₂ 28% EtOH
-53°C	=	Eutectic mixture EtOH-H ₂ O
-23°C	=	Liquid solid - CCl ₄

after reaching the bath temperature, ampoules are quickly transferred to -80°C bath, stored 24 hrs.

Selected Experimental Results:

† Waxman tube ampoules

I	Immersion bath	% recovery	Original motility
	-190° (liq N ₂)	9	30%
	-79°C	57	30%
	Cooled 1°/min.	25%	30%
II	-160°C	8%	25%
	-79°C	51%	25%
	-23°C	22%	25%

~~IV~~ Annular Ampoules:

Immersion bath	% Recovery	Original motility
-160°C	7%	47%
-80°C	18%	47%
-53°C	22%	47%
-23°C	13%	47%

Speed measurements, 25-24°C

No. moving sperm counted	Sample	Speed	mean dev.
260	II (control)	31 μ /sec	± 2
150	II (frozen, thawed)	29 μ /sec	± 2
365	III (control)	24 μ /sec	± 0.1
135	III (frozen, thawed)	21 μ /sec	± 2
275	I (control)	28 μ /sec	± 1.5
160	I (frozen, thawed)	25 μ /sec	± 2

Dear Szekely,

Marcus repeated the freezing expt. I give data below.

June 30, 1953

Assay before freezing 9.1×10^8
Mutant frequency before freezing
T4 $147/10^8$
T5 $7.1/10^8$
T6 $14.2/10^8$

	Control	Expt
Assay after freezing	—	11.2×10^8
Regrow 10cc sample in 1 liter broth. to about same density.		
Mutant frequency after regrow		
T4	$267/10^8$	$306/10^8$
T5	$19.3/10^8$	$16.7/10^8$
T6	$23.9/10^8$	$22.0/10^8$

Human sperm protocol was used.

I think the experiment was done as well as possible and I can think of no improvements

Caron